

Hardware & Software Requirements



Supported Operating System

Windows 10 Pro & Above

Minimum Hardware Requirements

CPU	Intel® Core™ i5-10500 Processor or similar
Video	Intel® UHD Graphics 630 or similar
Memory	8 GB RAM
Hard Drive	100 GB
Network	10/100/1000 Ethernet connection
Keyboard and mouse (Required)	USB Wired Keyboard and Wired Two-Button Mouse
Audio speakers (Optional)	External Stereo Speakers (2W)
Monitor (Required)	21" Color Monitor - 1920 x 1080 @ 60 Hz

Network

- Notepad ++ - Diagnostic Support Tool
- All CenTrak folders must be excluded from Anti-Virus scanning and live protection for optimum performance
- CenTrak needs read and write permissions to our working directory
 - C:\Program Files (x86)\Centrak
 - C:\ProgramData\Centrak\NotifierClient