

## Hardware & Software Requirements



### Supported Operating System

Windows 10 Pro

### Minimum Hardware Requirements

CPU	Intel® Core™ i5-10500 Processor or similar
Video	Intel® UHD Graphics 630 or similar
Memory	16 GB RAM
Hard Drive	500 GB
Network	10/100/1000 Ethernet connection
Keyboard and mouse	USB Wired Keyboard and Wired Two-Button Mouse
Audio speakers	External Stereo Speakers (2W)
Monitor	21" Color Monitor - 1920 x 1080 @ 60 Hz

### Network

- ASP .NET FRAMEWORK 4.0 or newer – supports the Operating System and Installers
- IIS 7 or newer - Internet Information Services (IIS) 7 Manager is an administration UI that provides end users and administrators with a way to remotely manage IIS 7.0 and IIS 7.5 servers.
- Windows Explorer 11 or Google Chrome browser – for our Security Solutions User Interface
- Adobe Reader - for reading Help files and saving Reports
- Adobe Flash\*\*
- Notepad ++ - Diagnostic Support Tool
- WinZip or similar compression software Diagnostic Support Tool
- All CenTrak folders must be excluded from Anti-Virus scanning and live protection for optimum performance

\*\*Needed for legacy CenTrak Systems only.