

## Hardware & Software Requirements



### Supported Operating System

Windows 10 Pro, Windows 11

### Minimum Hardware Requirements

CPU	Intel® Core™ i5-10500 Processor or similar
Video	Intel® UHD Graphics 630 or similar
Memory	16 GB RAM
Hard Drive	500 GB
Network	10/100/1000 Ethernet connection
Keyboard and mouse	USB Wired Keyboard and Wired Two-Button Mouse
Audio speakers	External Stereo Speakers (2W)
Monitor	21" Color Monitor - 1920 x 1080 @ 60 Hz

### Network

- Edge or Google Chrome browser – for our Safety Solutions User Interface
- Adobe Flash\*\*
- Notepad ++ - Diagnostic Support Tool
- All CenTrak folders must be excluded from Anti-Virus scanning and live protection for optimum performance
- Allow ICMP and port 80 and 443 between the Workstation IP and SS App Server IP

\*\*Needed for legacy CenTrak Systems only.

### Dedicated Workstation and Power Requirements:

- 24/7 operation
- No auto log out
- No sleep timer
- No hibernation timer
- No display sleep
- Auto login into the OS